



Item # _____

Job _____

MetroMax i™ I.V. TRANSPORT/STORAGE SYSTEM

The Perfect Solution for:

- Exchange System
- Par Level Storage
- Bulk Storage in Central Supply
- Pharmacy Admixture
- E.R., I.C.U., and O.R. Storage

- **Comprehensive System:** Six models of carts are constructed from a combination of state-of-the-art polymer components and epoxy-coated open-wire baskets and dividers.
- **Easily Accessible:** Baskets have semi-open front for easy access. Baskets fit level on frame for bulk storage or at a 10° angle for stock rotation. Dividers adjust horizontally within the basket to accept various sizes of I.V. bags.
- **Organized/Identifiable:** Dedicate an entire basket or a basket section by I.V. bag content or size. Continuous label holders identify baskets contents. Open-wire construction promotes quick visual inventory.
- **Corrosion Resistant:** The I.V. Transport/Storage System superstructure is rust and corrosion proof.
- **Strong and Durable:** Reinforced polymers give the cart's superstructure incredible strength, which enhances rigidity and eliminates swaying found in other I.V. storage and transport systems.
- **Easily Maneuverable:** 5" (127mm) casters, two with brakes, roll easily over thresholds and around corners.
- **Improved Infection Control:** MetroMax i™ helps improve infection control with a cleaner design that includes corrosion proof polymer materials, smooth rounded corners to allow for easy cleaning and *Microban® Antimicrobial Product Protection that keeps the cart cleaner between cleanings.



MXIV4



MXIV1



Metro incorporates several elements in its product design to support a facility's infection control processes:

- Microban® antimicrobial protection helps keep products "cleaner between cleanings"
- Advanced polymer and other proprietary finishes provide corrosion resistance
- Smooth rounded corners to allow for easier cleaning

Look for the "red check" symbol for this added protection



*Microban® antimicrobial product protection helps keep surfaces "cleaner between cleanings" by inhibiting the growth of bacteria, mold and mildew that cause odors and stains on the shelf surface.

*MICROBAN® and the MICROBAN® symbol are registered trademarks of the Microban Products Company, Huntersville, NC.




InterMetro Industries Corporation
North Washington Street
Wilkes-Barre, PA 18705
www.metro.com



MetroMax i™ I.V. Transport/Storage System

30.20


MetroMax i™ I.V. TRANSPORTATION/STORAGE SYSTEM



IBV1 (Dividers sold separately)



MXIV7



Baskets fit on frames level or at a 10° angle

Dimensions

Model No.	Description	Width		Length		Height		Total Baskets Per Cart	Casters
		(in.)	(mm)	(in.)	(mm)	(in.)	(mm)		
MXIV1	Open Cart	24	610	24	610	60	1524	4	*
MXIV2	Open Cart with Top Shelf	24	610	24	610	69	1753	4	*
MXIV4	Open Cart	24	610	42	1066	60	1525	8	†
MXIV5	Open Cart with Top Shelf	24	610	42	1066	69	1753	8	†
MXIV6	Open Cart	24	610	60	1524	60	1525	12	†
MXIV7	Open Cart with Top Shelf	24	610	60	1524	69	1753	12	†

NOTE: Each basket sold with cart includes two dividers.

* Four 5" (127mm) poly casters; two swivel, two with brake.

† 5" (127mm) poly casters on an aluminum dolly frame; two swivel, one swivel with brake, and one rigid.

Accessories

Model No.	Description	Size (Width x Depth x Height)	
		(in.)	(mm)
IBV1	I.V. Basket with Label Holder	18½ x 24½ x 10	470 x 622 x 254
IVBD	Divider	9	229

Epoxy-coated baskets and dividers have built-in Microban antimicrobial product protection.

All Metro Catalog Sheets are available on our Web Site: www.metro.com



InterMetro Industries Corporation
North Washington Street, Wilkes-Barre, PA 18705
Phone: 570-825-2741
Fax: 570-825-2852

L02-148e
Printed in U.S.A. Rev. 6/10
Information and specifications are subject to change without notice. Please confirm at time of order.

Copyright © 2010 InterMetro Industries Corp.

For Product Information:
U.S. and Canada: 1.800.992.1776

Latin America: 1.561.333.3824
Europe: +31.76.587.7550

Asia/Pacific: +65.6567.8003
Middle East/Africa: +971.4.811.8286